



Reagent set monochloramine, (indophenol method) HR TNT

Produktnr:

2805145

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.



Föreslagna ersättningar:

- 2802246

For determination of Monochloramine by the Indophenol Test 'N Tube[™] method. Hach Method 10172. Range: 0.1 - 10.0 mg/L as Cl₂. Sample Size: 2 mL. Up to 50 tests/reagent set.

Test 'N Tube products provide safe, convenient testing

Capped 16-mm vials provide a self-contained package for mixing and measurement.

Easy ordering

All necessary reagents and vials are contained in the package.

Specifikationer

Antal tester :	50
Fotnot :	+ Dependent on the number of performed reagent blanks.
Instrument:	All except DR 900, PC II
Mätområde:	0.1 - 10.0 mg/L
Metodnamn:	Indophenol
Metodnummer:	10172
Parameter:	Chloramine, Mono
Plattform:	Test 'N Tube