



APA6000 Alkalinity Reagent 1, Acid Titrant, 1 L

Produktnr:

2826153

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



The APA6000 Alkalinity Process analyser uses five reagents each of which is supplied in a 1 L bottle designed to provide reagents for 30 days of continuous use. APA Alkalinity Reagent 1 (2826153) is an acid titrant that reacts with the alkalinity in the sample. APA Alkalinity Reagent 2 (2696653) is an indicator used to determine the endpoints of the titration. APA Alkalinity Standards 1 and 2 (6001100) are used together to perform a multiple point calibration. The fifth reagent, APA Alkalinity Cleaning Solution (2697053), cleans the system during the Prime and Instrument Clean functions.

Specifikationer

Instrument:	APA6000
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Mätområde:	N/A
Parameter:	Alkalinity

Obligatoriska tillbehör

- Skyddshandskar storlek L, blå, nitril, puderfria, 50 styck (Item SM743L)