



Bottle, BOD, 300 mL, 24/pk, serialized: 1-24

Produktnr:	2898700
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included Serialized bottles have silk-screened, sequential numbers (1-24) and include a bar code.

Specifikationer

Body Material:	Glass
Kapacitet :	300 mL
Material:	Borosilicate Glass
Package size UOM:	24/pk, Serialized: 1 - 24
Volym/förpackningsstorlek:	24 /pk