



Bottle, storage, autoclavable, wide mouth, 125 mL pk/12

Produktnr:	2939103
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Polyethylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121°C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Specifikationer

Body Material:	Plastic
Kapacitet :	125 mL
Material:	PP
Package size UOM:	12pk
Volym/förpackningsstorlek:	12 /pk