



Bottle BOD, 300 mL, disposable

Produktnr:	2943101
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Single 300 mL disposable BOD bottle without stopper. Plastic bottle is EPA approved for BOD measurement.

Box of 100 bottles available under article number 2943100. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Specifikationer

Body Material:	Plastic
Kapacitet :	300 mL
Material:	PET
Volym/förpackningsstorlek:	1 each