



SPADNS 2 (arsenikfri) fluoridreagenslösning, 1 L

Produktnr:

2947553

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For determination of Fluoride by Hach SPADNS 2 Method 10225. The SPADNS2 reagent formulation removes interference from chlorine without arsenic. Range: 0.02 - 2.00 mg/L F⁻. EPA compliant. 1 Liter bottle.

Identical performance to traditional Hach SPADNS

Use the same Hach instruments, test procedures, and calibration curves

Specifikationer

Antal tester :	250
EPA compliant:	Ja
Mätområde:	0.02 - 2.00 mg/L F ⁻
Metodnamn:	SPADNS 2
Metodnummer:	10225
Parameter:	Fluoride
Plattform:	Lösning