



MultiTest kit case (145 x 206 x 66 mm), blue polypropylene

Produktnr:	4660100
SEK Pris:	Kontakta oss
Tillgänglig	

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Specifikationer

Dimension (H x B x D):	145 mm x 206 mm x 66 mm
Höjd :	5.75 in
Kapslingstyp:	B