



QL filters, pore size 8 - 12 μm , diameter 55 mm, circular, 100/ pk, qualitative

Produktnr: 50652

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specifikationer

Diameter:	55 mm
Material:	Cellulose
Volym/förpackningsstorlek:	100 /pk