



QL filters, pore size 8 - 12 µm, diameter 90 mm, circular, 100/ pk, qualitative

Produktnr:	50655
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specifikationer

Diameter:	90 mm
Material:	Cellulose
Volym/förpackningsstorlek:	100 /pk