



Cylinder, graduated, 5 mL $\pm$ 0.1 mL, 0.1 mL divisions (white markings)

Produktnr:

50837

SEK Pris:

Kontakta oss

Tillgänglig

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes $\geq$ 25 mL.

Specifikationer

Certifieringar :	Class B
Höjd :	125 mm
Kapacitet :	5 mL
Material:	Glass
Mätmetod:	To Deliver
Typ:	Single Metric
Volym/förpackningsstorlek:	1 each