



Cylinder, graduated, 10 mL +- 0.1 mL, 0.2 mL divisions (white markings)

Produktnr: 50838
SEK Pris: Kontakta oss
Inget leveransdatum anges

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specifikationer

Certifieringar :	Class B
Höjd :	135 mm
Kapacitet :	10 mL
Material:	Glass
Mätmetod:	To Deliver
Typ:	Single Metric
Volym/förpackningsstorlek:	1 each