



Cylinder, graduated, 50 mL +/- 0.4 mL, 1.0 mL divisions (white markings)

Produktnr:	50841
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specifikationer

Certifieringar :	Class B
Höjd :	190 mm
Kapacitet :	50 mL
Material:	Glass
Mätmetod:	To Contain
Typ:	Single Metric
Volym/förpackningsstorlek:	1 each