



**Be Right™**



## **B4M TOC Analyser 6 Month Service Kit**

|                   |                   |
|-------------------|-------------------|
| <b>Produktnr:</b> | <b>19-KIT-125</b> |
| SEK Pris:         | Kontakta oss      |
| Tillgänglig       |                   |

For BioTector B7000 MR analysers supplied from July 2014 to August 2015 with  
analyser serial number format BAC-4xxx-xxxx

---

### **Vad finns i förpackningen**

Includes: Omron Relays (3), EMPP Tubing, Y Tube fitting (2), SR25 pump tube rail (2), PTFE Diaphragm, CO2 Filter, Ozone Destructor catalyst, Filter Mats for vent (2), Lens Tissue for CO2 Analyser, Viton o-ring (2)