



Recharge kit of 4 pre-filled cartridges with premounted 2956A membranes for A1100 oxygen sensors with 33051-xT cap only

Produktnr:

2956A-AT05

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Includes o-rings and cleaning tools. Special low concentration electrolyte.

Specifikationer

Instrument:

Orbisphere A1100-S00T O2 sensor, Orbisphere GA2x00-S00T O2 sensor

Lagringsförhållanden:

Store in a dry place at room temperature

Parameter:

O2