



Sension+ 5209 pH-vätskefylld kombinationselektrod för mikroprovsvolymer, 6 mm diameter, skruvlock S7

Produktnr: LZW5209.97.0002
SEK Pris: Kontakta oss
Tillgänglig

pH liquid electrode (refillable) with ceramic diaphragm. Glass body. Range 0-14 pH, 0 to 80 °C. Screw cap S7. 6 mm diameter, 8 mm minimum immersion depth.

This pH electrode is suitable for Micro Sample volumes.

For Micro Sample volumes

High stability reference electrode due to sleeved cartridge

Encapsulated reference system

Specifikationer

Accuracy:	0.002 pH with pH31 Meter
Dimension (D x L):	8 mm x 100 mm
Elektrolyt:	LZW9500.99 KCl 3M
Garanti :	12 månader
Givar koppling:	Ceramic
Givartyp:	Glass
Kontakt:	Screw cap S7
Längd:	100 mm
Material:	Sensor Body: Glass
Mätområde:	0 - 14 pH
Parameter:	pH
Temperaturområde :	0 - 80 °C

Obligatoriska tillbehör

- AS7-kabel / 1M / BNC för instrument med BNC-kontakt (Item LZW9055.99)