



Be Right™



Orthophosphate Standard Solution, 0 mg/L PO₄-P, 500 mL

| | |
|-------------------|----------------|
| Produktnr: | LCW1051 |
| SEK Pris: | Kontakta oss |
| Tillgänglig | |

Standard Solution for use with the NP6000sc Phosphate Analyser.

Specifikationer

| | |
|----------------------------|---------------------------|
| Instrument: | NP6000sc |
| Koncentration: | 0 mg/L PO ₄ -P |
| Parameter: | Orthophosphate |
| Volym/förpackningsstorlek: | 500 mL |