



**Be Right™**



## **Orthophosphate Standard Solution, 5 mg/L PO<sub>4</sub>-P, 500 mL**

<b>Produktnr:</b>	<b>LCW1054</b>
SEK Pris:	Kontakta oss
Tillgänglig	

Standard Solution for use with the NP6000sc Phosphate Analyser.

---

### **Specifikationer**

Instrument:	NP6000sc
Koncentration:	5 mg/L PO <sub>4</sub> -P
Parameter:	Orthophosphate
Volym/förpackningsstorlek:	500 mL