



Be Right™



Chlorine, free; colour disc test kit, model CN-66F (pk/100)

Produktnr:

223102-EU

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Measuring range 0.1 - 3.4 mg/L Cl₂

Specifikationer

Antal tester :	100
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kapslingstyp:	B
Mätområde:	Single Range
Metodnamn:	Color disc/DPD
Modell:	CN-66F
Parameter:	Chlorine Free
Plattform:	Color Disc
Smallest increment:	0.2
Standardmätområde:	0 - 3.4 mg/L Cl ₂
Test-kit modell:	CN-66F ²