



Test kit, nitrogen, nitrate color disc, model NI-14

Produktnr:

1416100-EU

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO₃ and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	50
Kapslingstyp:	B
Mätområde:	Dual Range
Metodnamn:	Color disc/ Cadmium reduction
Minsta ökningssteg:	0,05
Modell:	NI-14
Parameter:	Nitrogen, Nitrate - As NO ₃ ⁻ -N
Plattform:	Color Disc
Smallest increment:	0,05 0,5
Standardmätområde:	0 - 1 mg/L NO ₃ ⁻ -N; 0 - 10 mg/L NO ₃ ⁻ -N
Test-kit modell:	NI-14**