



Nitrate test kit, model NI-11

Produktnr:

146803-EU

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Farligt



One of Hach's most popular test kits, it uses the simple, low-cost color disc method.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	100
Kapslingstyp:	B
Mätområde:	0 - 40 mg/L NO ₃ ⁻ -N
Metodnamn:	Color disc/ Cadmium reduction
Modell:	NI-11
Parameter:	Nitrogen, Nitrate - As NO ₃ ⁻ -N
Plattform:	Color Disc
Smallest increment:	1.0
Standardmätområde:	0 - 40 mg/L NO ₃ ⁻ -N
Test-kit modell:	NI-11**